



# Baby630

PROFESSIONAL ELECTRONIC BABY SCALE



Folding baby tray  
(W)280x(L)360x(H)175



BABY630: CARRY CASE

| Model   | Code  | Capacity | Division | Dimension           |
|---------|-------|----------|----------|---------------------|
| Baby630 | 00488 | 7.5/15kg | 5/10g    | Total               |
|         |       |          |          | W)560x(L)380x(H)160 |
|         |       |          |          | Tray                |
|         |       |          |          | (W)560x(L)240       |

## FUNCTIONS

Weight, Tare, Hold weight, Weight percentile, Reset-to-zero, AutoOff programmable.

## PRODUCT DETAILS

|                                  |   |
|----------------------------------|---|
| <b>Certifications:</b>           | Directive NAWI Approved class III. ECC directive MDD Risk class Im, with measurement function |
| <b>Structure:</b>                | ABS display and platform; Folding baby tray ABS.  |
| <b>Display:</b>                  | Large LCD 25mm 5 digit easy to reading.   |
| <b>Power:</b>                    | External adapter or 6 alkaline batteries code AA 1,5V.  |
| <b>Weight:</b>                   | Approx 3.5kg  |
| <b>Accessories included:</b>     | Level Bubble, Adjustable Feet, Adapter, 6 alkaline batteries code AA 1,5V                     |
| <b>Accessories not included:</b> | Carry case cod.00157C.  |
| <b>Packaging:</b>                | 1pcs dimension (w)450x(L)200x(h)300mm, weight approx 4kg                                      |
| <b>Warranty:</b>                 | 2 years.  |