

MIMLED 600

Compact LED Surgery lamp. MIMLED 600 offers 60.000 lux at 1 meter.

HIGH COLOUR RENDERING

The high colour rendering $Ra > 94\%$, describes the exceptional quality of white light, allows tissues to be distinguished with utmost precision.

LIGHT INTENSITY DIMMING

Adjustable light intensity integrated in the handgrip and on the head which can be regulated using the panel control.

ENERGY SAVING LIGHT

Low power consumption.

25W MIMLED 600

AVAILABLE IN 3 OPTIONS

- Trolley Stand.
- Wall version.
- Ceiling version (up to 3.10 meter height).



MIMLED 600

| Technical Details | |
|----------------------|----------------------|
| Lux at 1 meter | 60.000 lux |
| Lamp lifetime | 50.000 hours |
| Colour temperature | 4.500 °K |
| Colour rendering | ≥ 94% |
| Light field diameter | 195 mm |
| Power supply | 100/230 V - 50/60 Hz |
| Consumption | 25 W |



Adjustable light intensity integrated in the grip



Ceiling version
Cod. ML600C



Wall version
Cod. ML600W



Trolley Stand version
Cod. ML600FL